PLASMA HIRUDIN STANDARD HIGH

SC020L

For *in vitro* use only

For research use only

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Last revision/ dernière révision : 12/01/2006

Intended use

Lyophilised plasmas, at **0** and **5 µg/mL** hirudin concentration, for establishing the calibration curve for the clotting assay of hirudin in human plasmas (Hemoclot Hirudin, #CK002K).

Reagents:

Each kit contains:

- Standard 0: Std 0: 3 vials of 1 mL of normal plasma, lyophilised, without hirudin (**0 µg/mL**)
- Standard 5: Std 5: 3 vials of 1 mL of normal plasma, lyophilised, containing **5 µg/mL** of hirudin.

Reagent Preparation and stability:

In their original package, and before any use, when stored at 2-8°C, the reagents are stable until the expiration date printed on the kit.

Excipients:

Glycine, stabilizers.

Preparation:

- **Std 0: 0 µg/mL** plasma:
 - restore the vial with 1mL distilled water
 - let for 15 min at room temperature
 - mix gently until complete dissolution of the content.
- **Std 5: 5 µg/mL** plasma:
 - restore the vial with 1mL distilled water
 - let for 15 min at room temperature
 - mix gently until complete dissolution of the content.

These plasmas are stable for:

- 24 hours at room temperature (18-25 °C)
- 48 hours at 2-8 °C
- 2 months, frozen in their original vial, or in a plastic tube (before use, thaw in a water bath at 37 °C, for at least 15 min).

Nota: Plasmas used for hirudin standards (0 and 5 µg/mL) were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Protocol for the calibration curve:

The Hemoclot hirudin assay uses a 20 fold (**1:20**) plasma dilution for the high range determination. Using the plasmas at **0 µg/mL** and **5 µg/mL** hirudin concentration, prepare the following calibration curve, using physiological saline (8,5 g/l Sodium Chloride) or Owren Koller type buffer as diluent. The calibration standards, already diluted 20 fold (**1:20**) are obtained as follows:

Plasma Hirudin concentration (µg/mL)	0	1.25	2.50	3.75	5
Plasma 0 µg/mL (µL)	100	75	50	25	0
Plasma 5 µg/mL (µL)	0	25	50	75	100
Physiological Saline or Owren Koller Buffer (µL)	1900	1900	1900	1900	1900

Diluted calibrators are stable for 6 hours at room temperature.

Complementary information

Plasma hirudin calibrators are standardised according to hirudin concentration. The hirudin protein activity can present variations from lot to lot, according to its specific activity (usually in the range 14,000 to 16,000 ATU*/mg)

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