BIOPHEN CS-51(09)

Factor IXa Chromogenic substrate

Vial of 25 mg/flacon de 25 mg

Ref. 229051

STORE AT
2-8 °C

For research use only
Not for use in diagnostic
procedures.

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CONSERVER A
2-8 °C

Origin:

Synthetic Chromogenic substrate, specific for Factor IXa, highly purified and stabilized.

Specificity:

Recommended substrate for the chromogenic measurement of factor IXa.

Presentation:

Vial containing 25 mg of BIOPHEN CS-51(09), lyophilised in presence of Mannitol as a bulking agent.

Reconstitution:

According to the research protocol used, the BIOPHEN CS-51(09) chromogenic substrate can be restored with variable volumes of distilled water; for example 5 mL can be used for a substrate concentration of 5 mg/mL, or 20 ml for a substrate concentration of 1.25 mg/mL, etc....

Peptide sequence:

CH₃SO₂-(D)-CHG-Gly-Arg-pNa, AcOH

Purity grade:

> 95%.

Molecular weight:

628.7.

Free pNA content (measured on substrate at 2.5g/L)

OD_{405nm} ≤ 0.30

Use:

All research studies and protocols where a source of chromogenic substrate for Factor IXa is required. Factor IXa chromogenic activity is dramatically enhanced by 33% Ethylene Glycol. **This kit is for research use only and should not be used for patient diagnosis or treatment.**

Respective reactivities:

FIXa	aPC	FXa	Plasmin	Thrombin
100	8	40	70	50

Assay conditions must be established for making the substrate totally specific for factor IXa.

Stability:

Restored :

- 7 days at room temperature
- 1 month at 2-8 °C
- **Do not freeze.**

Lyophilized : Until the expiration date printed on the vial.

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