Purified Human Fibrinogen, Plasminogen Free Fibrinogène humain purifié sans Plasminogène

1 vial (PP002A) / 1 flacon (PP002A)
Set of 6 vials (PP002K) / Coffret de 6 flacons (PP002K)
REF PP002A / PP002K

**FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC PROCEDURES.**

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English, last revision 02-2018

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INTENDED USE:

All research studies and protocols where a source of highly purified fibrinogen is necessary.
This kit is for research use only and should not be used for patient diagnosis or treatment.

REAGENTS:

[Pn] Highly purified fibrinogen extracted from citrated human plasma.
Plasminogen is specifically removed by affinity chromatography.
Vial containing at least 4.5 mg of highly purified human fibrinogen, plasminogen free.
1 vial of 1 mL (PP002A).
6 vials of 1 mL (PP002K).

WARNING AND PRECAUTIONS:

- The plasma used for fibrinogen purification was tested with registered methods and found negative for HIV antibodies, HBs Ag and HCV antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.
- Waste should be disposed of in accordance with applicable local regulations.
- Handle the reagents with care to avoid contamination during use. If possible, avoid reagent evaporation during use by limiting the liquid-air exchange surface.
- Aging studies show that the reagents can be shipped at room temperature without degradation.
- To preserve reagent stability, seal the vials after use with their respective caps.
- This device is intended for professional use in the laboratory.
- For *in vitro* use.

REAGENT PREPARATION AND STABILITY:

Reconstitute the contents of each vial with exactly **1 mL distilled water**, in order to obtain a 5 mg/mL fibrinogen solution, shake vigorously until fully dissolved. This solution can be diluted to the desired concentration when required. Allow to stabilize for 15 min. at room temperature (18-25°C), shaking occasionally.
Homogenize the reagent prior to use.
Reagent stability after reconstitution, excluding any contamination or evaporation, and stored in the original vial, is of:

- **96 hours** at 2-8°C.
- **48 hours** at room temperature (18-25°C).

STORAGE CONDITIONS:

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

LIMITATIONS:

- Any reagent presenting an unusual appearance or showing signs of contamination must be rejected.
- **Excipients:** - Glycine, trisodium citrate and HEPES buffer.
- 6 mg sodium chloride/vial.
- **Purity:** One major band of 340,000 daltons on SDS-PAGE. This fibrinogen has a clottability \geq 98%.
- **Activity:** Tested for the clotting activity induced by thrombin. When the fibrinogen, at 5 mg/mL (100 μ L), is clotted by thrombin, at 5 NIH/mL (100 μ L), the clotting time must be below 15 sec. Absence of plasminogen is tested with a clot lysis method, in presence of tPA.
The exact activity is reported on the analysis certificate for each substrate.

The results obtained should be for research purposes only and not used for patient diagnosis or treatment.

SYMBOLS:

Symbols used and signs listed in the ISO 15223-1 standard, see Symbol definitions document.

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