Purified Human Fibrinogen, Plasminogen Free

REF PP002A, Pn 1 x 1 mL

REF PP002K, Pn 6 x 1 mL

Purified Human Fibrinogen, Plasminogen Free
FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC PROCEDURES.

English, last revision 07-2020

INTENDED USE:

All research studies and protocols where a source of highly purified fibrinogen is necessary.

This kit is for research use only and should not be used for patient diagnosis or treatment.

REAGENTS:

Pn Highly purified fibrinogen extracted from citrated human plasma.

Plasminogen is specifically removed by affinity chromatography.

Vial containing at least 4.5 mg of highly purified human fibrinogen, plasminogen free.

REF PP002A : 1 vial of 1 mL

REF PP002K : 6 vials of 1 mL

- **Excipients:** - Glycine, trisodium citrate and HEPES buffer.
- 6 mg sodium chloride/vial.
- **Purity:** One major band of 340,000 daltons on SDS-PAGE. This fibrinogen has a clotability ≥ 98%.
- **Activity:** Tested for the clotting activity induced by thrombin. When the fibrinogen, at 5 mg/mL (100 µL), is clotted by thrombin, at 5 NIH/mL (100 µL), the clotting time must be below 15 sec. Absence of plasminogen is tested with a clot lysis method, in presence of tPA

WARNING AND PRECAUTIONS:

- Some reagents provided in these kits contain materials of human origin. Whenever human plasma is required for the preparation of these reagents, approved methods are used to test the plasma for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.
- Waste should be disposed of in accordance with applicable local regulations.
- Aging studies show that the reagents can be shipped at room temperature without degradation.
- This device of *in vitro* use is intended for professional use in the laboratory.

REAGENT PREPARATION:

Gently remove the freeze-drying stopper, to avoid any product loss when opening the vial.

Pn Reconstitute the contents of each vial with exactly **1 mL of distilled water**, in order to obtain a 5 mg/mL fibrinogen solution, shake vigorously until fully dissolved. This solution can be diluted to the desired concentration when required. Allow to stabilize for 15 min. at room temperature (18–25°C), shaking occasionally. Homogenize the reagent prior to use.

STORAGE AND STABILITY:

Unopened reagents should be stored at 2–8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

Pn Reagent stability after reconstitution, free from any contamination or evaporation, and stored closed, is of:

- **96 hours** at 2–8°C.
- **48 hours** at room temperature (18–25°C).

LIMITATIONS:

- Any reagent presenting an unusual appearance or showing signs of contamination must be rejected.

The results obtained should be for research purposes only and not used for patient diagnosis or treatment.

SYMBOLS:

Symbols used and signs listed in the ISO 15223-1 standard, see Symbol definitions document.

Changes compared to the previous version.

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