CRYO*check*[™] **IVD**

LA CHECK[™]

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CRYO*check*[™]

dRVVT Screening Reagent

LA CHECK[™]

Intended Use

CRYO*check* LA Check is a dilute Russell’s Viper Venom Time (dRVVT) reagent intended to screen for the presence of lupus anticoagulants (LA) in citrated human plasma.

Summary and Principle

LA are autoantibodies of the IgG and IgM types that are specifically directed against negatively charged phospholipids, such as phosphatidylinositols and phosphatidylserines, or complexes of phospholipids with either β₂-glycoprotein-1 or clotting factors such as prothrombin¹. They occur in various clinical conditions, especially autoimmune diseases². LA have traditionally been detected using phospholipid sensitive *in vitro* clotting tests, such as the activated partial thromboplastin time (APTT), kaolin clotting time (KCT), and dRVVT³. The dRVVT was introduced in 1986 and showed improved sensitivity to LA over the APTT partially due to a reduced phospholipid concentration⁴.

LA are usually indicated by a prolonged clotting time result that is not corrected by mixing patient plasma with normal plasma. The correction of a prolonged result by the addition of phospholipids to the plasma is a more specific characteristic of LA⁴.

LA prolong phospholipid-sensitive clotting tests; however, they are paradoxically associated with thrombotic problems⁶. LA are a common cause of unexplained prolonged APTTs and need to be carefully distinguished from idiopathic antibodies against factor VIII associated with bleeding.

LA are now considered to be a significant risk factor in patients with otherwise unexplained thrombosis and are often present in women who have recurrent fetal loss⁶. They are also associated with a variety of hemostatic problems such as thrombocytopenia and neurological disorders⁷.

Russell’s viper venom directly activates factor X, bypassing factor VII of the extrinsic pathway and the contact and antihemophilic factors of the intrinsic pathway. Therefore, dRVVT tests are more specific for LA than APTTs as they are not affected by contact factor abnormalities or by factor VIII deficiencies or antibodies³.

Reagents

CRYO*check* LA Check contains Russell’s viper venom, phospholipids, antiheparin agents, calcium, buffers, stabilizers, sodium azide, and green dye (<0.001%).

 *Sodium azide may react with lead and copper plumbing to form highly explosive metal compounds. Ensure proper disposal of reagent according to federal, state, and local regulations.*

Storage and Handling

When stored at -40 to -80°C, CRYO*check* LA Check is stable to the end of the month indicated on the product packaging.

Thaw each vial at 37°C (± 1°C) in a waterbath. **The use of a dry bath or heating block for thawing is not recommended.** Thawing times are important and should be strictly adhered to. The use of a timer is recommended. Refer to the Thawing Table for recommended thawing duration based on aliquot size. Allow thawed reagent to acclimate to room temperature (18 to 25°C) and invert gently prior to use.

Thawing Table	
Aliquot Size	37°C (± 1°C) Waterbath
1.0 mL	4 minutes

CRYO*check* LA Check may be used for up to 48 hours after thawing, if capped in the original vial and maintained at 2 to 8°C. Allow refrigerated reagent to acclimate to room temperature (18 to 25°C) and invert gently prior to use. **Thawed material may be refrozen once, and stored at -20°C for up to one month.**

Availability

Product	Catalog #	Format
LA Check	CHK-10	25 vials x 1.0 mL

Instruments

Each lab should prepare the local instrument in accordance with the manufacturer’s instructions for use.

Procedure

Materials Provided

- CRYO*check* LA Check

Materials Required but not Provided

- Waterbath capable of maintaining 37°C (± 1°C)
- 12 mm x 75 mm glass test tubes
- Quality control material (e.g. CRYO*check* Lupus Positive Control, CRYO*check* Weak Lupus Positive Control)
- Plastic disposable pipettes
- Volumetric pipette
- Timer
- Stopwatch

Specimen Collection and Preparation

Patient samples should be collected into 105 - 109 mmol/L sodium citrate dihydrate anticoagulant (3.2%) in a ratio of 9 parts blood to 1 part anticoagulant. Patient plasma is derived by centrifugation at 1500 x g for 15 minutes in order to achieve platelet-poor plasma (<10,000 platelets/μL), and should be tested within four hours of collection when maintained at 2 to 4°C in accordance with CLSI guidelines⁸. If samples are to be frozen before testing, plasmas should be centrifuged a second time, and stored at -20°C or below.

Manual Method – Tilt Tube

- In a 37°C (± 1°C) waterbath, prewarm a slight excess of CRYO*check* LA Check allowing 200 μL per test.
- Dispense 200 μL of test plasma into a test tube and warm for one minute at 37°C (± 1°C).
- Add 200 μL of prewarmed CRYO*check* LA Check to the plasma and simultaneously initiate the clot timer. Record clotting times in seconds.
- Repeat for duplicate test values and report the average of these as the result.

Automated Methods

Reagent preparation instructions and instrument settings for a variety of analyzers are available upon request from Precision BioLogic.

Quality Control

Each laboratory should establish its own quality control (QC) ranges using acceptable statistical methods. These QC ranges may then be used to monitor and validate the integrity of the test system⁹. For all coagulation tests, the laboratory must include at least two levels of control for every eight hours of operation and any time a change in reagents occurs¹⁰.

Commercial lyophilized quality control plasmas containing unspecified levels of citrate and platelets are not recommended as they may give erroneous results^{11,12}.

Results

If the CRYO*check* LA Check clotting time is within the laboratory-established normal reference range, further testing for LA may not be necessary. If the CRYO*check* LA Check clotting time result is prolonged (i.e. more than 2 standard deviations (SD) longer than the mean of the laboratory-established normal reference range), the result should be investigated further. A prolonged result may be indicative of:

- the presence of LA
- factor II, V, or X deficiency
- oral anticoagulant therapy (OAT)

If the CRYO*check* LA Check result is prolonged, CRYO*check* LA Sure[™] confirmatory test (Precision BioLogic Catalog No. SUR-10) should be performed. Evidence of inhibitory activity can be shown using an additional mixing step of patient plasma and pooled normal plasma. Failure of the mixture to correct the original prolonged clotting time is evidence of a possible inhibitor. This step may be incorporated into the initial CRYO*check* LA Check screening procedure⁵.

In accordance with the SSC Subcommittee for the Standardization of LA guidelines⁵, results should be compared to other established LA tests performed on the same sample, since no single assay can guarantee, with certainty, that LA is present or absent.

Limitations of the Procedure

Patients with deficiencies of factors II, V, or X or patients on anti-vitamin K therapy may exhibit prolonged CRYO*check* LA Check times. Although normal plasma mixing studies may correct for these deficiencies, published data has demonstrated that further dilution of weak non-specific inhibitors such as a lupus anticoagulant can occur and produce false-negative results in dRVVT procedures¹³. CRYO*check* LA Check results on plasmas subjected to mixing studies should be interpreted with care.

CRYO*check* LA Check is unaffected by heparin levels up to 1.0 unit/mL. Plasmas containing heparin levels greater than 1.0 unit/mL may give false-positive results and should not be tested with this reagent.

Plasma samples with visible hemolysis should not be used due to possible clotting factor activation and endpoint measurement interference⁸. Icteric or lipemic samples may also interfere with endpoint determination on some optical instruments⁸.

Expected Values

In a study of 20 healthy males and females using a Diagnostica Stago ST4[®] analyzer, a CRYO*check* LA Check normal reference range (2 SD confidence interval) of 29.1 - 42.7 seconds was established. These values should be used as a guide only. Each laboratory should establish their own normal reference range.

Performance Characteristics

In precision studies over 48 hours at 2 to 8 C with CRYO*check* Lupus Positive Control plasma on a Diagnostica Stago ST4[®] analyzer, CRYO*check* LA Check exhibited an overall coefficient of variation (CV) of 3.26%.

An R²=0.969 was derived in a correlation study using Gradipore LA Screen[™] dRVVT screening test involving OAT patient plasmas (n=15), known LA positive samples (n=12), and plasmas with depleted levels of factors II (n=3), V (n=3), and X (n=3).

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Bibliography / Bibliographie

- Arnout J. Antiphospholipid syndrome: Diagnostic aspects of lupus anticoagulants. Thromb Haemost 2001; 86:83-91.
- Feinstein DI. Lupus anticoagulant, thrombosis and fetal loss. N Eng J Med 1985; 313:1248-1350.
- Triplett DA, Brandt JT. Laboratory identification of the lupus anticoagulant. Br J Haematol 1991; 73:139-142.
- Thiagarajan P, Pengo V, Shapiro SS. The use of the dilute Russell’s Viper Venom Time for the diagnosis of lupus anticoagulants. Blood 1986; 68:869-874.
- Brandt JT, Triplett DA, Alving B, Scharrer I. Criteria for the diagnosis of lupus anticoagulants: an update. On behalf of the Subcommittee on Lupus Anticoagulant/Antiphospholipid Antibody of the Scientific and Standardisation Committee of the ISTH. Thromb Haemost 1995 Oct; 74(4):1185-90
- Love PE, Santoro SA. Antiphospholipid antibodies: Anticardiolipin and the lupus anticoagulant in systemic lupus erythematosis (SLE) and in non-SLE disorders. Ann Intern Med 1990; 112:682-698.
- Harris EN, Asherson RA, Hughes GRV. Antiphospholipid antibodies – autoantibodies with a difference. Ann Revs Med 1998; 39:261-271.
- Collection, transport, and processing of blood specimens for coagulation testing and general performance of coagulation assays, NCCLS, H21-A3. 1998.
- Cembrowski GS, Carey RN. Laboratory quality management. Chicago: ASCP Press; 1989.
- CLIA 1988 – Code of Federal Regulations, 42CFR493.1253, 1998.
- Adcock DM, Kressin DC, Marlar RA. Effect of 3.2% vs. 3.8% sodium citrate concentration on routine coagulation testing. Am J Clin Pathol 1997; 107:105-10.
- Hirst CF, Poller L. The cause of turbidity in lyophilized plasmas and its effects on coagulation tests. J Clin Pathol 1992; 45:701-3.
- Thom J, Ivey L, Eikelboom J. Normal plasma mixing studies in the laboratory diagnosis of lupus anticoagulant. J Thromb Haemost 2003; 1: 2689-91

	European Authorized Representative (Regulatory affairs only) Emergo Europe Molenstraat 15, 2513 BH The Hague, The Netherlands Tel: (+31) 70 345 8570 Fax: (+31) 70 346 7299
<i>Symbols used / Symboles utilisés</i>	09.60.00017 Rev. 13 February / février 2011
<i>In vitro</i> diagnostic medical device	Batch code
	Use by
	Temperature limitation
	Warning
	Manufacturer
	Authorized representative
	Désignation du lot
<i>in vitro</i> diagnostic	Date de péremption
	Températures limites de conservation
	Attention
	Fabricant
	Mandataire

<i>In vitro</i> diagnostic medical device	Batch code	Use by	Temperature limitation	Warning	Manufacturer	Authorized representative
Dispositif médical <i>in vitro</i> diagnostic	Désignation du lot	Date de péremption	Températures limites de conservation	Attention	Fabricant	Mandataire