



# Human Plasminogen Matched Pair Antibodies for EIA (5 x 96 Tests)

**REF** PG-EIA

For Research Use Only  
Not for Use in Diagnostic Procedures  
For in vitro use only

**Store at -10 to -20° C**

## INTENDED USE:

Human Plasminogen Matched Pair Antibodies for EIA are intended for use with in-house enzyme-linked immunosorbent assays for measuring human Plasminogen in plasma, or in any biological fluid where human Plasminogen can be present. The results obtained should be for research purposes only and not used for patient diagnosis or treatment.

## SUMMARY:

Plasminogen is synthesized in the liver and circulates in plasma at a concentration of ~190 µg/mL (~2.3 µM). Plasminogen is a single-chain glycoprotein of ~88 kDa that consists of a catalytic domain followed by five kringle structures. Within these kringle structures are four low-affinity lysine binding sites and one high-affinity lysine binding site. It is through these lysine binding sites that plasminogen binds to fibrin and to α2-Antiplasmin. Native Plasminogen (Glu-Plasminogen) exists in two variants that differ in their extent of glycosylation, and each variant has up to six isoelectric forms with respect to sialic acid content, for a total of 12 molecular forms.

Activation of Glu-Plasminogen by the Plasminogen activators Urokinase (uPA), or tissue Plasminogen Activator (tPA) occurs by cleavage after residue Arg560 to produce the two-chain active serine protease Plasmin. In a positive feedback reaction, the Plasmin generated cleaves a ~8 kDa peptide from Glu-Plasminogen, producing Lys77-Plasminogen which has a higher affinity for Fibrin and when bound is a preferred substrate for Plasminogen activators such as Urokinase. Additional activators of Plasminogen include Kallikrein and activated Factor XII.

## ASSAY PRINCIPLE:

The diluted plasma sample or biological fluid is introduced into one of the microwells of a micro ELISA plate which has been pre-coated with anti-human Plasminogen antibody. When present in the added material, Plasminogen binds to the anti-human polyclonal antibody. Following a washing step, the remaining bound antibodies are revealed with a detection antibody, anti-human Plasminogen peroxidase conjugate, which reacts specifically with human Plasminogen. Following another washing step, the peroxidase substrate, o-Phenylenediamine (OPD) in presence of hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>), is introduced and a yellow color develops. The color turns orange when the reaction is stopped with sulfuric acid. The color developed is directly proportional to the amount of Plasminogen present in the tested sample.

## REAGENTS:

### Required Materials provided (enough for 5x96 Tests):

- **C:** Capture Antibody (PG-EIA-C). 1 vial of 0.5 mL polyclonal affinity-purified antibody specific for Plasminogen. For coating plates. Yellow cap.
- **D:** Detecting Antibody (PG-EIA-D). 1 vial of 0.5 mL polyclonal antibody specific for human Plasminogen, coupled to peroxidase. For detecting captured Plasminogen. Red cap.

**Note:** Antibodies are provided in a glycerol solution (50% v/v) and should be stored at **-10 to -20°C**. Vials should be tightly capped. Do not store in frost-free freezers

Antibodies can be centrifuged briefly in a micro-centrifuge to gather residual reagent from the cap and walls of the tube.

In their original packaging, before use, when stored at -10 to -20°C, the unopened antibodies are stable until the expiration date printed on the vial.

### Required Materials not provided:

Optimum performance can be obtained when the following solutions and assay conditions are used.

- **Micro ELISA plates** with hydrophilic surface designed for high binding of IgG. For example, 96-well Immulon 4-HBX.
- **Coating Solution** (50mM Carbonate). Dissolve 1.59 g of Na<sub>2</sub>Cl<sub>3</sub> and 2.93 g of NaHCO<sub>3</sub> in distilled water to a final volume of 1 L and adjust pH to 9.6. Store at 2-8°C for 1 month.
- **Phosphate-Buffered Saline (PBS)** [For preparation of wash and blocking solutions.] Dissolve 8.0 g NaCl, 1.15 g Na<sub>2</sub>HPO<sub>4</sub>, 0.2 g KH<sub>2</sub>PO<sub>4</sub> and 0.2 g KCl in distilled water to a final volume of 1 L and adjust pH to 7.4. Store up to 1 month at 2-8°C, discard if there is evidence of microbial growth.
- **Wash Solution:** (PBS/Tween-20 0.1% v/v). Add 1.0 mL of Tween-20 to 1 L of PBS and adjust pH to 7.4. Store at 2-8°C up to 1 week.
- **Blocking Solution:** (PBS/BSA 1% w/v). Dissolve 2.5 g of Bovine Serum Albumin in 200 mL of PBS and adjust pH to 7.4; add PBS to final volume of 250 mL. Aliquot and store frozen at -20°C.
- **Sample Diluent:** (HEPES/BSA/Tween-20). Dissolve 5.95 g HEPES (free acid), 1.46 g NaCl and 2.5 g Bovine Serum Albumin in 200 mL distilled H<sub>2</sub>O; add 0.25 mL of Tween-20 and adjust pH to 7.2 with NaOH; add distilled water to final volume of 250 mL. Aliquot and store frozen at -20°C for up to 6 months.
- **Substrate Solution:** (Citrate-Phosphate buffer). Dissolve 2.6 g Citric Acid and 6.9 g Na<sub>2</sub>HPO<sub>4</sub> in 450 mL distilled H<sub>2</sub>O and adjust pH to 5.0 with Phosphoric acid or NaOH; add distilled water to final volume of 500 mL. Store at 2-8°C up to 1 month.
- **OPD Substrate**(o-phenylenediamine.2HCl) ☒ Toxic! 5 mg tablets: e.g. Sigma #P-6912. Prepare immediately before use. Dissolve 5 mg OPD in 12 mL Substrate Solution and then add 12 µL 30% H<sub>2</sub>O<sub>2</sub>. Do not store.
- **Stop Solution:** (2.5M H<sub>2</sub>SO<sub>4</sub>) ☒ Corrosive! Generates heat on dilution! Handle with great care. Avoid any skin and eye contact. Wear protective glasses and gloves when handling. Carefully add 13.9 mL 18M H<sub>2</sub>SO<sub>4</sub> to 86 mL distilled H<sub>2</sub>O. Store at room temperature for up to 1 month.
- **Reference standards** for Plasminogen which have the same matrix and anticoagulant as the samples to be tested
- Micro ELISA plate washing equipment and shaker.
- Plate reader with a wavelength set up at 490 nm.

**PROCEDURE:**

**1. Coat ELISA plate:** Dilute the Capture Antibody with Coating Solution 1/100 (use polypropylene tube) and immediately add 100  $\mu$ L to every well in the plate. Incubate for 2 hours at 22°C or overnight at 2-8°C.

**2. Blocking:** Empty contents of plate and add 150  $\mu$ L of Blocking Solution to every well and incubate for 60 minutes at 22°C. This step blocks any remaining binding sites on the plastic wells. Wash plate 3X with Wash Solution.

**3. Samples:** Dilute Plasminogen Reference standard with Sample Diluent 1/10,000 (100%) then serially dilute by halves down to 1/320,000 (3.13%). Dilute sample plasmas or biological fluid with Sample Diluent 1/20,000, 1/40,000 and 1/80,000. Apply 100  $\mu$ L per well and incubate plate at 22°C for 60 minutes. Wash plate 3X with Wash Solution. (Plasma samples should not be applied at dilutions lower than 1/20, as falsely high readings may result.)

**4. Detecting Antibody:** Dilute the Detecting Antibody with Sample Diluent 1/100 and apply 100  $\mu$ L to each well. Incubate plate at 22°C for 60 minutes. Wash plate 3X with Wash Solution.

**5. OPD Substrate:** Apply 100  $\mu$ L of freshly prepared OPD substrate to each well. Allow color to develop for 5-10 minutes then stop color reaction with the addition of 50  $\mu$ L per well of Stop Solution. Read the plate at a wavelength of 490 nm. (Optimal color development time is the time required to obtain A490  $\geq$  1.000 for the 100% reference point, not to exceed 20 minutes.)

**Additional Notes:**

- Do not shake plate during incubation.
- Do not allow the wells to become dry. Keep plate covered or in a humid chamber during incubations.
- Rheumatoid factor in samples may bind to the capture and/or detecting antibodies and cause interference in the ELISA assay.

**6. Calibration Curve:** On bi-logarithmic graph paper, plot the known Plasminogen concentrations on abscissa and the corresponding absorbance (A490) on ordinates in order to establish the calibration curve.

**RESULTS:**

From the constructed calibration curve, directly determine the Plasminogen concentration and multiply by the appropriate dilution factor.

---

CoaChrom Diagnostica GmbH | Your Coagulation Experts  
Hauptstrasse 5  
2344 Maria Enzersdorf, Austria  
Phone: +43-1-236 222 1 | Fax: +43-1-236 222 111  
Phone: 0800-246 633 0 | Fax: 0800 - 246 633 3 (Toll free)  
info@coachrom.com | www.coachrom.com