**Tris-NaCl-EDTA Buffer**  
**pH 7.85**  
**1x20 mL / Set of 4x20 mL**  
**REF. AR032A / AR032K**

**Tris-NaCl-EDTA Puffer**  
**pH 7,85**  
**1x20 mL / Set mit 4x20 mL**  
**Art.Nr. AR032A / AR032K**

For *in-vitro* diagnostic use

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**INTENDED USE:**

Special dilution buffer reducing Heparin Interferences for anti-Xa chromogenic assays, especially using BIOPHEN™ DiXal kit (221030) or BIOPHEN™ Heparin LRT kits (221011/221013/221015).

**REAGENTS:**

**[Bf]** **Tris-NaCl-EDTA buffer, pH 7.85.** Ready to use. Contains 1% PEG and a heparin neutralizing substance.

**1 vial of 20 mL (AR032A)**

**4 vials of 20 mL (AR032K)**

Reagent contains small amounts of sodium azide (0.9 g/L), see WARNINGS AND PRECAUTIONS.

**WARNINGS AND PRECAUTIONS:**

- In contact with lead or copper pipes, sodium azide can generate explosive compounds.
- Waste should be disposed of in accordance with applicable local regulations.
- Handle the reagents with care to avoid contamination during use. If possible, avoid reagent evaporation during use by limiting the liquid-air exchange surface. Evaporation reduces the reagent's stability in the analyzer.
- Aging studies show that the reagents can be shipped at room temperature without degradation.
- To preserve reagent stability, seal the vials after use with their respective caps.
- For *in vitro* diagnostic use.

**REAGENTS PREPARATION AND STABILITY:**

**[Bf]** Ready to use. Allow to stabilize for 30 minutes at room temperature (18-25°C) before use.

Homogenize the reagent prior to use, while avoiding formation of foam.

Reagent stability after opening, excluding any contamination or evaporation, and stored in the original vial, is of:

- **2 months** at 2-8°C
- **7 days** at room temperature (18-25°C)

**STORAGE CONDITIONS:**

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

**LIMITATIONS:**

- Any reagent presenting an unusual appearance or showing signs of contamination must be rejected.
- Refer to the specific associated kit inserts for instructions concerning the tests to be performed.

**SYMBOLS:**

Symbols used and signs listed in the ISO 15223-1 standard, see Symbol definitions document.

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