# Special Calcium Chloride

## Heparin neutralizing

**REF** AR036K **S-Ca** 6 vials x 10 mL / 6 flacons x 10 mL

Special Calcium Chloride

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English, revision: 11-2025

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### INTENDED USE:

Special Calcium Chloride solution with heparin neutralization in activated Partial Thromboplastin Time (aPTT) assays using an automated method. This device of *in vitro* diagnostic use is intended for professional use in the laboratory.

### REAGENTS:

**S-Ca** **Special Calcium Chloride** Solution at approximately 0.025 mol/L, ready to use. Contains a heparin neutralizing substance and preservative.

The product is classified as non-hazardous and is not subject to labelling according to EC Regulation No. 1272/2008 [CLP].

### WARNING AND PRECAUTIONS:

- Waste should be disposed of in accordance with applicable local regulations.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### REAGENT PREPARATION:

**S-Ca** Reagent is ready to use, homogenize and load it directly on the analyzer following Application Guide instruction.

### STORAGE AND STABILITY:

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

**S-Ca** Reagent stability after opening, free from any contamination or evaporation, and stored closed, is of:

- 60 days at 2-8°C.
- Stability on board of the analyzer: see the specific Application Guide.

### REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED:

- For use with aPTT reagents manufactured by Hyphen Biomed Only on CS-series, ACL-TOP® family, CN-series.

### LIMITATIONS:

- Any reagent presenting no limpid appearance or showing signs of contamination must be rejected.
- Refer to the specific associated kit inserts for instructions concerning the test to be performed.

### RESULTS:

- The aPTT is not affected by Heparin (Unfractionated heparin and low molecular weight heparin) up to 1 IU/mL.
- Lot to lot variability measured on 3 lots is: %CV ≤10%.

eIFU is available on the HYPHEN BioMed website.

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