

## Certificate of Analysis

### Plasminogen (Pg) Deficient Human Plasma

<b>Article #:</b>	PG-MP
<b>Lot / Exp.:</b>	SAMPLE
<b>Preparation:</b>	Pooled normal citrated human plasma depleted of Pg using antibodies directed to plasminogen immobilized on agarose beads. Plasma contains 20 mM HEPES.
<b>Storage:</b>	Plasma is shipped frozen and should be stored below -60°C.
<b>Presentation:</b>	Frozen Pg deficient plasma.
<b>Specifications:</b>	Assay results:

Assay	Lot# xxx
PT (HemosIL PT-Fib)	14.8 sec
APTT (HemosIL APTT-SP)	33.0 sec
Fibrinogen (clottable)	3.37 g/L
Plasminogen (antigen) by ELISA	<0.01 U/ml

Although this material was prepared with plasma collected from donors screened for CJD and which was tested at source and found negative for HBsAg, RPR and antibodies to HIV and HCV and non-reactive for HIV-1 rNA and HCV rNA by FDA approved tests, it should be handled by personnel trained in proper procedures for handling potential viral contaminants.